Strep A Rapid Test Device (Throat Swab)

English: Package Insert

Page 1

A rapid test for the qualitative detection of Strep A antigen in throat swab specimens.

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Strep A Rapid Test Device
(Throat Swab)
Package Insert
English

INTENDED USE

The Strep A Rapid Test Device (Throat Swab) is a rapid chromatographic immunoassay for the qualitative detection of Strep A antigen from throat swab specimens to aid in the diagnosis of Group A Streptococcal infection.

A rapid test for the qualitative detection of Strep A antigen in throat swab specimens. For professional *in vitro* diagnostic use only.

SUMMARY

Streptococcus pyogenes is non-motile gram-positive cocci, which contains the Lancefield Group A antigen that can cause serious infections such as pharyngitis, respiratory infection, impetigo, endocarditis, meningitis, puerperal sepsis, and arthritis.¹ Left untreated, these infections can lead to serious complications, including rheumatic fever and peritonsillar abscess.² Traditional identification procedures for Group A Streptococci infection involve the isolation and identification of viable organisms using techniques that require 24 to 48 hours or longer.^{3,4}

The Strep A Rapid Test Device (Throat Swab) is a rapid test to qualitatively detect the presence of Strep A antigen in throat swab specimens, providing results within 5 minutes. The test utilizes antibodies specific for whole cell Lancefield Group A *Streptococcus* to selectively detect Strep A antigen in a throat swab specimen.

PRINCIPLE

The Strep A Rapid Test Device (Throat Swab) is a qualitative, lateral flow immunoassay for the detection of Strep A carbohydrate antigen in a throat swab. In this test, antibody specific to Strep A carbohydrate antigen is coated on the test line region of the test. During testing, the extracted throat swab specimen reacts with an antibody to Strep A that is coated onto particles. The mixture migrates up the membrane to react with the antibody to Strep A on the membrane and generate a color line in the test line region. The presence of this color line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region. If the control line does not appear, the test result is not valid.

REAGENTS

The test contains Strep A antibody coated particles and Strep A antibody coated on the membrane.

PRECAUTIONS

- For professional *in vitro* diagnostic use only. Do not use after expiration date.

- Each test device is for single use only.
- Standard personal hygiene measures should be taken in the case of ingestion or direct eye contact with the reagent.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed
- Avoid touching finger directly to the specimen pad or result window of the test as this can cause incorrect results.
- The used test should be discarded according to local regulations.
- Humidity and temperature can adversely affect results.
- The test must remain in the sealed pouch until use.
- Do not use test if pouch is damaged.
- Need to read the results in the required time.
- Reagent B contains an acidic solution. If the solution contacts the skin or eye, flush with large volumes of water.
- The positive and negative controls contain sodium azide (NaN_3) as a preservative.
- Do not mix reagents from different lots.
- Do not interchange reagent bottle caps.
- Do not interchange external control solution bottle caps.
- Proper storage condition is critical to product performance.
- Once the foil pouch is opened, use the product as soon as possible.
- Maybe the bystander is contaminated with the biological material and sample in the testing process, so stay away from the bystander.
- Test the product in the prescribed environmental condition.
- Use the correct specimen type.
- Do not recycle the device.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

- Only use reagents and sterile swabs provided in the kit.
- Collect the throat swab specimen with the sterile swab that is provided in the kit. The sterile swabs provided with this test must be used for specimen collection. Transport swabs containing modified Stuart's or Amies medium can also be used with this product. Swab the posterior pharynx, tonsils and other inflamed areas. Avoid touching the tongue, cheeks and teeth with the swab.⁵
- Testing should be performed immediately after the specimens have been collected. Swab specimens may be stored in a clean, dry plastic tube for up to 8 hours at room temperature

(15-30°C) or 72 hours at 2-8°C.

- If a culture is desired, lightly roll the swab tip onto a Group A selective (GAS) blood agar plate before using the swab in the Strep A Rapid Test Device (Throat Swab).
- Specimens should not be frozen and thawed repeatedly.

MATERIALS

Materials Provided

• Test devices • Sterile swabs • Dropper tips	• Workstation • Package insert • Test tubes
• Strep A Reagent A (13.8% Sodium Nitrite)	 Warning Harmful if swallowed. Wash thoroughly after handling. Do not eat, drink or smoke when using this product. IF SWALLOWED: Call a POISON Center or doctor/physician if you feel unwell. Rinse mouth. Dispose of contents/container in according with local regulation.
• Strep A Reagent B (1.2% Acetic Acid) • Strep A Positive control (Non-viable Strep A; 0.09% NaN₃) • Strep A Negative control (Non-viable Strep C; 0.09% NaN₃)	Safety data sheet available for professional users on request

Materials Required But Not Provided

- Timer


DIRECTIONS FOR USE

Allow the test, reagents, throat swab specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Remove the test device from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.
2. Hold the Reagent A bottle vertically and **add 4 full drops** (approximately 240 µL) of Reagent A to an extraction test tube. Reagent A is red in color. Hold the Reagent B bottle vertically and **add 4 full drops** (approximately 160 µL) to the tube. Reagent B is colorless. Mix the solution by gently swirling the extraction test tube. The addition of Reagent B to Reagent A changes the color of the solution from red to pale yellow. See illustration 1 & 2.
3. **Immediately add the throat swab to the extraction test tube of yellow solution.** Agitate the swab **10 times** in the tube. Leave the swab in the tube for **1 minute**. Then press the swab against the side of the tube and squeeze the bottom of the tube as the swab is withdrawn. Discard the swab. See illustration 3.
4. Fit the dropper tip on top of the extraction test tube. Place the test device on a clean and

level surface. **Add 3 full drops of solution** (approx. 100 μ L) to the specimen well (S) and then start the timer. See illustration 4.

5. Wait for the colored line(s) to appear. **Read the result at 5 minutes.** Do not read the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE: * **Two lines appear.** One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T). A positive result indicates that Strep A was detected in the specimen.

***NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of Strep A present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

NEGATIVE: **One colored line appears in the control line region (C).** No line appears in the test line region (T). A negative result indicates that Strep A antigen is not present in the specimen, or is present below the detectable level of the test. The patient's specimen should be cultured to confirm the absence of Strep A infection. If clinical symptoms are not consistent with results, obtain another specimen for culture.

INVALID: **No line appears in the control line region (C).** If this occurs, read the directions again and repeat the test with a new test. If the result is still invalid, stop using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal Quality Control

Internal procedural controls are included in the test. A color line appearing in the control region (C) is an internal positive procedural control.

External Quality Control

It is recommended that a positive and negative external control be run every 25 tests, and as deemed necessary by internal laboratory procedures. External positive and negative controls are supplied in the kit. Alternatively, other Group A and non-Group A *Streptococcus* reference strains may be used as external controls. Some commercial controls may contain interfering preservatives; therefore, other commercial controls are not recommended.

Procedure for External Quality Control Testing

1. Add 4 full drops of Reagent A and 4 full drops of Reagent B into an extraction test tube. Tap the bottom of the tube gently to mix the liquid.
2. Add 1 full drop of positive or negative control solution into the tube, holding the bottle upright.
3. Place a clean swab into this extraction tube and agitate the swab in the solution by rotating it at least 10 times. Leave the swab in the extraction tube for 1 minute. Then express the liquid from the swab head by rolling the swab against the inside of the extraction tube and squeezing the extraction tube as the swab is withdrawn. Discard the swab.
4. Continue with Step 4 of Directions For Use.

If the controls do not yield the expected results, do not use the test results. Repeat the test or contact your distributor.

LIMITATIONS

1. The Strep A Rapid Test Device (Throat Swab) is for *in vitro* diagnostic use only. The test should be used for the detection of Strep A antigen in throat swab specimens only. Neither the quantitative value nor the rate of increase in Strep A antigen concentration can be determined by this qualitative test.
2. This test will only indicate the presence of Strep A antigen in the specimen from both viable and non-viable Group A *Streptococcus* bacteria.
3. A negative result should be confirmed by culture. A negative result may be obtained if the concentration of the Strep A antigen present in the throat swab is not adequate or is below the detectable level of the test.
4. Except Puritan Polyester Sterile Swab, Changfeng Polyester Sterile Swab, BBL Amies Liquid Swab and BBL Stuart Liquid Swab, other swabs have not been validated to use with this test.
5. Excess blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks, and teeth⁵ and any bleeding areas of the mouth with the swab when collecting specimens.
6. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.

EXPECTED VALUES

Approximately 15% of pharyngitis in children ages 3 months to 5 years is caused by Group A beta-hemolytic *Streptococcus*.⁶ In school-aged children and adults, the incidence of Strep throat infection is about 40%.⁷ This disease usually occurs in the winter and early spring in temperate climates.³

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

Using three medical centers for evaluation, a total of 492 throat swabs were collected from patients exhibiting symptoms of pharyngitis. Each swab was rolled onto a sheep blood agar plate, and then tested by the Strep A Rapid Test Device (Throat Swab). The plates were further streaked for isolation, and then incubated at 37°C with 5-10% CO₂ and a Bacitracin disk for 18-24 hours. The negative culture plates were incubated for an additional 18-24 hours. Possible GAS colonies were subcultured and confirmed with a commercially available latex agglutination grouping kit.

Of the 492 total specimens, 384 were confirmed to be negative and 108 were confirmed to be positive by culture. During this study, two Strep F specimens yielded positive results with the Test. One of these specimens was re-cultured, then re-tested and yielded a negative result. Three additional different Strep F strains were cultured and tested for cross-reactivity and also yielded negative results.

Method		Culture		Total Results
Strep A Rapid Test	Results	Positive	Negative	
	Positive	102	7	109
	Negative	6	377	383
Total Results		108	384	492

Relative Sensitivity: 94.4% (88.3%-97.9%)*

Relative Specificity: 98.2% (96.3%-99.3%)*

Accuracy: 97.4% (95.5%-98.6%)*

* 95% Confidence Intervals

Precision

Intra-Assay

Within-run precision has been determined by using three levels of specimen: a low positive, a medium positive and a high positive value specimens. The results were correctly identified 100% of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the low positive specimen. Two different lots of the Strep A Rapid Test Device (Throat Swab) have been tested using a low positive, a medium positive and a high positive value specimens. The specimen is correctly identified 100% of the time.

Cross-Reactivity

The following organisms were tested at 1.0×10^7 organisms per test and were all found to be negative when tested with the Strep A Rapid Test Device (Throat Swab). No mucoid-

producing strains were tested.

Group B Streptococcus
Group F Streptococcus
Streptococcus pneumoniae
Streptococcus mutans
Staphylococcus aureus
Corynebacterium diphtheria
Candida albicans
Pseudomonas aeruginosa

Neisseria meningitidis
Neisseria sicca
Branhamella catarrhalis
Group C Streptococcus
Group G Streptococcus
Streptococcus sanguis
Enterococcus faecalis
Staphylococcus epidermidis

Serratia marcescens
Klebsiella pneumoniae
Bordetella pertussis
Neisseria gonorrhea
Neisseria subflava
Hemophilus influenza

Interfering Substances

The Strep A Rapid Test Device (Throat Swab) has been tested for specimens containing Blood, 1% Cepacol, 1% Cheryl halls, 1% Dimetapp, 1% Listerine, 1% menthol Halls, 1mg Mucin, 1% Robitussin, 1% Scope and 1% Vicks. Results indicate that no interference was observed in these specimens.

POL Studies

Three physicians' offices were used to conduct an evaluation of the Strep A Rapid Test Device (Throat Swab). Personnel with various educational backgrounds performed the testing. Each physician's office tested a randomly coded panel of samples consisting of negative (20), low positive (20), and medium positive (20) for three days. The results obtained had a 100% (95% CI: 98.3% - 100%) correlation with the expected results.

Index of Symbols

	Consult instructions for use		Contains sufficient for <n> tests		Authorized representative in the European Community
	In vitro diagnostic medical device		Use-by date		Do not reuse
	Store between 2-30°C		Batch code		Catalogue number
	Warning		Manufacturer		Reagent A
	Reagent B		Positive Control		Negative Control

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BIBLIOGRAPHY/BIBLIOGRAPHIE/BIBLIOGRAFIA/BIBLIOGRAFIA

1. Murray, P.R., et al. Manual of Clinical Microbiology, 6th Edition, ASM Press, Washington D.C., 1995, p. 299-307.
2. Webb, KH. *Does Culture Confirmation of High-sensitivity Rapid Streptococcal Tests Make Sense? A Medical Decision Analysis*. Pediatrics (Feb 1998), 101:2, 2.
3. Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. *Diagnosis and Management of Group A Streptococcal Pharyngitis*. Clinical Infectious Diseases (1997), 25: 574-83.
4. Needham CA, McPherson KA, Webb KH. *Streptococcal Pharyngitis: Impact of a High-sensitivity Antigen Test on Physician Outcome*. Journal of Clinical Microbiology (Dec 1998), 36: 3468-3473.
5. Shea, Y.R., *Specimen Collection and Transport*, Clinical Microbiology Procedures Handbook, Isenberg, H.D., American Society of Microbiology, Washington D.C., 1.1.1-1.1.30, 1992.
6. Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. *Group A beta-hemolytic streptococcal pharyngitis in preschool children aged 3 months to 5 years*. Clinical Pediatrics (June 1999), 38: 357-360.
7. Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA. *Group A Streptococcal Pharyngitis in Adults 30 to 65 years of age*. Southern Medical Journal (May 1999), 491-492.

Abon Biopharm (Hangzhou) Co., Ltd.
#198 12th Street East,
Hangzhou Economic & Technological
Development Area,
Hangzhou, 310018, P.R.China
www.abon.com.cn

EC REP

MT Promedt Consulting GmbH
Altenhofstrasse 80
66386 St. Ingbert/Germany
Tel: +49 (0) 6894 581020

Effective date:2020-06-24

Number: 1156136102